

兔抗 UAP1 多克隆抗体

中文名称:	兔抗 UAP1 多克隆抗体
英文名称:	Anti-UAP1 rabbit polyclonal antibody
别名:	AGX; AGX1; AGX2; SPAG2
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	UAP1
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	UAP1 belongs to the DPGP type 1 family. It converts UDP and GlcNAc-1-P into UDP-GlcNAc, and UDP and GalNAc-1-P into UDP-GalNAc. Isoform AGX1 has 2 to 3 times higher activity towards GalNAc-1-P, while isoform AGX2 has 8 times more activity towards GlcNAc-1-P.
Applications:	ELISA, WB, IHC
Name of antibody:	UAP1
Immunogen:	Fusion protein of human UAP1
Full name:	UDP-N-acetylglucosamine pyrophosphorylase 1
Synonyms:	AGX; AGX1; AGX2; SPAG2
SwissProt:	Q16222
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100
WB Predicted band size:	59 kDa
WB Positive control:	Human testis tissue lysate
WB Recommended dilution:	200-1000



