

## 兔抗 UBAC1 多克隆抗体

中文名称:	兔抗 UBAC1 多克隆抗体
英文名称:	Anti-UBAC1 rabbit polyclonal antibody
别名:	GBDR1; UBADC1
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
宿主:	Rabbit
抗原:	UBAC1
反应种属:	Human Mouse Rat
标记物:	Unconjugate
克隆类型:	Unconjugate

### 技术规格

<b>Background:</b>	Non-catalytic subunit of the KPC complex that acts as E3 ubiquitin-protein ligase. Required for poly-ubiquitination and proteasome-mediated degradation of CDKN1B during G1 phase of the cell cycle.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	UBAC1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human UBAC1.
<b>Full name:</b>	UBA domain containing 1
<b>Synonyms :</b>	GBDR1; UBADC1
<b>SwissProt:</b>	Q9BSL1
<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	293 cells and Huvec cells lysates

WB Recommended dilution: 500-3000



