

兔抗 UBAP2L 多克隆抗体

中文名称:	兔抗 UBAP2L 多克隆抗体
英文名称:	Anti-UBAP2L rabbit polyclonal antibody
别名:	NICE-4
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
宿主:	Rabbit
抗原:	UBAP2L
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	a ubiquitous protein containing an UBA-domain, which are present in multiple proteins of the ubiquitination pathway. Encoded by a gene of the human epidermal differentiation complex (EDC). May be involved in epidermal differentiation. Three isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	UBAP2L
Immunogen:	Synthesized peptide derived from internal of human UBAP2L.
Full name:	ubiquitin associated protein 2-like
Synonyms :	NICE-4
SwissProt:	Q14157

WB Predicted band size:	115 kDa
WB Positive control:	HT29 cells and 293 cells lysates
WB Recommended dilution:	500-3000

