

兔抗 UBXN2A 多克隆抗体

中文名称：兔抗 UBXN2A 多克隆抗体

英文名称：Anti-UBXN2A rabbit polyclonal antibody

别名：UBXD4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UBXN2A

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein is located in both the ER and cis-Golgi compartments from biochemical and immunofluorescence studies.UBXD4 as one of the cytosolic proteins that interact directly with the $\alpha 3$ and $\alpha 4$ nAChR subunits.Overexpression of UBXD4 in differentiated PC12 cells (dPC12) increased nAChR cell surface expression, especially that of the $\alpha 3\beta 2$ subunit.
Applications:	ELISA, WB
Name of antibody:	UBXN2A
Immunogen:	Synthetic peptide of human UBXN2A
Full name:	UBX domain-containing protein 2A
Synonyms :	UBXD4
SwissProt:	P68543

ELISA Recommended dilution:	1000-2000
WB Predicted band size:	29 kDa
WB Positive control:	Mouse liver tissue
WB Recommended dilution:	200-1000

