

## 兔抗 UBXN10 多克隆抗体

中文名称: 兔抗 UBXN10 多克隆抗体

英文名称: Anti-UBXN10 rabbit polyclonal antibody

别名: UBX domain protein 10; UBXD3

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: UBXN10

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Depletion of UBXN10 in zebrafish embryos causes defects in left-right asymmetry, which depends on functional cilia.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBXN10
<b>Immunogen:</b>	Synthetic peptide of human UBXN10
<b>Full name:</b>	UBX domain protein 10
<b>Synonyms:</b>	UBXD3
<b>SwissProt:</b>	Q96LJ8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	20-100

