

兔抗 UCK2 多克隆抗体

中文名称: 兔抗 UCK2 多克隆抗体

英文名称: Anti-UCK2 rabbit polyclonal antibody

别名: UK; UMPK; TSA903

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: UCK2

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a pyrimidine ribonucleoside kinase. The encoded protein (EC 2.7.1.48) catalyzes phosphorylation of uridine and cytidine to uridine monophosphate (UMP) and cytidine monophosphate (CMP), respectively.
Applications:	ELISA, IHC
Name of antibody:	UCK2
Immunogen:	Synthetic peptide of human UCK2
Full name:	uridine-cytidine kinase 2
Synonyms:	UK; UMPK; TSA903
SwissProt:	Q9BZX2
ELISA Recommended diluti	5000-10000

on:	
IHC positive control:	Human tonsil
IHC Recommend dilution:	20-100

