

兔抗 UGT3A1 多克隆抗体

中文名称: 兔抗 UGT3A1 多克隆抗体

英文名称: Anti-UGT3A1 rabbit polyclonal antibody

相关类别: 一抗

抗 原: UGT3A1

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	UDP-glucuronosyltransferases catalyze phase II biotransformation reactions in which lipophilic substrates are conjugated with glucuronic acid to increase water solubility and enhance excretion. They are of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.
Applications:	ELISA, WB, IHC
Name of antibody:	UGT3A1
Immunogen:	Fusion protein of human UGT3A1
Full name:	UDP glycosyltransferase family 3 member A1
SwissProt:	Q6NUS8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human gastric cancer

IHC Recommend dilution:	100-200
WB Predicted band size:	59 kDa
WB Positive control:	Human placenta tissue, Human between peritoneal stromal sarcoma tissue lysates
WB Recommended dilution:	500-2000

