

兔抗 UNC5A 多克隆抗体

中文名称: 兔抗 UNC5A 多克隆抗体

英文名称: Anti-UNC5A rabbit polyclonal antibody

级别: UNC5H1

储存: 冷冻 (-20℃)

抗原: UNC5A

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	UNC5A belongs to a family of netrin-1 (MIM 601614) receptors thought to mediate the chemorepulsive effect of netrin-1 on specific axons. For more information on UNC5 proteins, see UNC5C (MIM 603610).
Applications:	ELISA, WB, IHC
Name of antibody:	UNC5A
Immunogen:	Synthetic peptide of human UNC5A
Full name:	unc-5 netrin receptor A
Synonyms:	UNC5H1
SwissProt:	Q6ZN44
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and human tonsil

IHC Recommend dilution:	25-100
WB Predicted band size:	93 kDa
WB Positive control:	Human testis tissue lysate
WB Recommended dilution:	200-1000





