

兔抗 USE1 多克隆抗体

中文名称: 兔抗 USE1 多克隆抗体

英文名称: Anti-USE1 rabbit polyclonal antibody

别名: D12; P31; SLT1; MDS032

抗原: USE1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclona

技术规格

Background:	SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER.
Applications:	ELISA, WB
Name of antibody:	USE1
Immunogen:	Full length fusion protein
Full name:	unconventional SNARE in the ER 1
Synonyms:	D12; P31; SLT1; MDS032
SwissProt:	Q9NZ43
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	29 kDa
WB Positive control:	Human fetal liver tissue

WB Recommended dilution:

200-1000

