

兔抗 USP15 多克隆抗体

中文名称: 兔抗 USP15 多克隆抗体

英文名称: Anti-USP15 rabbit polyclonal antibody

别名: UNPH4; UNPH-2

储存: 冷冻 (-20℃)

抗原: USP15

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a member of the ubiquitin specific protease (USP) family of deubiquitinating enzymes. USP enzymes play critical roles in ubiquitin-dependent processes through polyubiquitin chain disassembly and hydrolysis of ubiquitin-substrate bonds. The encoded protein associates with the COP9 signalosome, and also plays a role in transforming growth factor beta signalling through deubiquitination of receptor-activated SMAD transcription factors. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 2.

Applications:	ELISA, IHC
Name of antibody:	USP15
Immunogen:	Fusion protein of human USP15
Full name:	ubiquitin specific peptidase 15
Synonyms :	UNPH4; UNPH-2
SwissProt:	Q9Y4E8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

