

兔抗 USP45 多克隆抗体

中文名称: 兔抗 USP45 多克隆抗体

英文名称: Anti-USP45 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃) 避光

宿 主: Rabbit

抗 原: USP45

反应种属: Human

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Belongs to the peptidase C19 family. 3 isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	USP45
Immunogen:	Synthesized peptide derived from internal of human USP45.
Full name:	ubiquitin specific peptidase 45
Synonyms :	0
SwissProt:	Q70EL2
WB Predicted band size:	92 kDa
WB Positive control:	A549 cells lysate
WB Recommended	500-3000

dilution:

