

兔抗 VGLL3 多克隆抗体

中文名称: 兔抗 VGLL3 多克隆抗体

英文名称: Anti-VGLL3 rabbit polyclonal antibody

别名: vestigial like family member 3; VGL3; VGL-3

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: VGLL3

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May act as a specific coactivator for the mammalian TEFs.
Applications:	ELISA, IHC
Name of antibody:	VGLL3
Immunogen:	Synthetic peptide of human VGLL3
Full name:	vestigial like family member 3
Synonyms:	VGL3; VGL-3
SwissProt:	A8MV65
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human gastric cancer
IHC Recommend dilution:	20-100

kDa

130—

95—

72—

55—

36—

28—

17—

