

兔抗 VPS51 多克隆抗体

中文名称: 兔抗 VPS51 多克隆抗体

英文名称: Anti-VPS51 rabbit polyclonal antibody

别名: FFR; ANG2; ANG3; C11orf2; C11orf3

相关类别: 一抗

抗原: VPS51

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the vacuolar protein sorting-associated protein 51 family. The encoded protein is a component of the Golgi-associated retrograde protein complex which acts as a tethering factor for carriers in retrograde transport from the early and late endosomes to the trans-Golgi network. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	VPS51

Immunogen:	Synthetic peptide of human VPS51
Full name:	VPS51 GARP complex subunit
Synonyms:	FFR; ANG2; ANG3; C11orf2; C11orf3
SwissProt:	Q9UID3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	40-200



