

兔抗 VSIG8 多克隆抗体

中文名称: 兔抗 VSIG8 多克隆抗体

英文名称: Anti-VSIG8 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: VSIG8

反应种属: Human, Mouse

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	VSIG8 contains 2 Ig-like V-type (immunoglobulin-like) domains. VSIG8 is single-pass type I membrane protein. The function of the VSIG8 protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	VSIG8
Immunogen:	Fusion protein of human VSIG8
Full name:	V-set and immunoglobulin domain containing 8
SwissProt:	P0DPA2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	25-100

