

兔抗 WDR45 多克隆抗体

中文名称: 兔抗 WDR45 多克隆抗体

英文名称: Anti-WDR45 rabbit polyclonal antibody

别名: JM5; NBIA4; NBIA5; WDRX1; WIPI4; WIPI-4

储存: 冷冻 (-20℃)

抗原: WDR45

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartate (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene has a pseudogene at chromosome 4q31.3. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene, but the biological validity and full-length nature of some variants

	ts have not been determined.
Applications:	ELISA, IHC
Name of antibody:	WDR45
Immunogen:	Fusion protein of human WDR45
Full name:	WD repeat domain 45
Synonyms :	JM5; NBIA4; NBIA5; WDRX1; WIPI4; WIPI-4
SwissProt:	Q9Y484
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

