

兔抗 WDR61 多克隆抗体

中文名称: 兔抗 WDR61 多克隆抗体

英文名称: Anti-WDR61 rabbit polyclonal antibody

别名: SKI8; REC14

储存: 冷冻 (-20℃)

抗原: WDR61

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

克隆类型: rabbit polyclonal

标记物: Unconjugate

技术规格

Background:	WDR61 is a subunit of the human PAF and SKI complexes, which function in transcriptional regulation and are involved in events downstream of RNA synthesis, such as RNA surveillance (Zhu et al., 2005 [PubMed 16024656]).
Applications:	ELISA, IHC
Name of antibody:	WDR61
Immunogen:	Full length fusion protein
Full name:	WD repeat domain 61
Synonyms :	SKI8; REC14
SwissProt:	Q9GZS3
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human esophagus cancer

IHC Recommend dilution:

20-100

