

## 兔抗 WFDC1 多克隆抗体

中文名称: 兔抗 WFDC1 多克隆抗体

英文名称: Anti-WFDC1 rabbit polyclonal antibody

别名: PS20

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: WFDC1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a member of the WAP-type four disulfide core domain family. The WAP-type four-disulfide core domain contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many family members. This gene is mapped to chromosome 16q24, an area of frequent loss of heterozygosity in cancers, including prostate, breast and hepatocellular cancers and Wilms' tumor. This gene is downregulated in many cancer types and may be involved in the inhibition of cell proliferation. The encoded protein may also play a role in the susceptibility of certain CD4 memory

	T cells to human immunodeficiency virus infection. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	WFDC1
<b>Immunogen:</b>	Synthetic peptide of humanWFDC1
<b>Full name:</b>	WAP four-disulfide core domain 1
<b>Synonyms :</b>	PS20
<b>SwissProt:</b>	Q9HC57
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	10-50

