

兔抗 WFDC5 多克隆抗体

中文名称: 兔抗 WFDC5 多克隆抗体

英文名称: Anti-WFDC5 rabbit polyclonal antibody

别名: PRG5; WAP1; dJ211D12.5

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: WFDC5

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the WAP-type four-disulfide core (WFDC) domain family. Most WFDC proteins contain only one WFDC domain, and this encoded protein contains two WFDC domains. The WFDC domain, or WAP signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor. Most WFDC gene members are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene belongs to the centromeric cluster.
Applications:	ELISA, IHC
Name of antibody:	WFDC5

Immunogen:	Fusion protein of human WFDC5
Full name:	WAP four-disulfide core domain 5
Synonyms :	PRG5; WAP1; dJ211D12.5
SwissProt:	Q8TCV5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human breast cancer
IHC Recommend dilution:	25-100



