

兔抗 XIRP1 多克隆抗体

中文名称：兔抗 XIRP1 多克隆抗体

英文名称：Anti-XIRP1 rabbit polyclonal antibody

别名：Xin; CMYA1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：XIRP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Protects actin filaments from depolymerization.
Applications:	ELISA, IHC
Name of antibody:	XIRP1
Immunogen:	Synthetic peptide of human XIRP1
Full name:	xin actin binding repeat containing 1
Synonyms:	Xin; CMYA1
SwissProt:	Q702N8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

