

兔抗 XKR7 多克隆抗体

中文名称: 兔抗 XKR7 多克隆抗体

英文名称: Anti-XKR7 rabbit polyclonal antibody

别名: C20orf159; dJ310O13.4

抗原: XKR7

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

XKR7 (XK-related protein 7) is a 579 amino acid multi-pass membrane protein that likely is a component of the XK/Kell complex of the Kell blood group system. The gene encoding XKR7 maps to human chromosome 20, which comprises approximately 2% of the human genome. Chromosome 20 contains nearly 63 million bases that encode over 600 genes, some of which are associated with Creutzfeldt-Jakob disease, amyotrophic lateral sclerosis, spinal muscular atrophy, ring chromosome 20 epilepsy syndrome and Alagille syndrome.

Applications:	ELISA, IHC
Name of antibody:	XKR7
Immunogen:	Synthetic peptide of human XKR7
Full name:	XK, Kell blood group complex subunit-related family, member 7
Synonyms :	C20orf159; dj310O13.4
SwissProt:	Q5GH72
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colorectal cancer and Human tonsil
IHC Recommend dilution:	25-100

