

兔抗 YY1AP1 多克隆抗体

中文名称: 兔抗 YY1AP1 多克隆抗体

英文名称: Anti-YY1AP1 rabbit polyclonal antibody

别名: YAP; HCCA1; HCCA2; YY1AP

储存: 冷冻 (-20℃)

抗原: YY1AP1

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	YY1AP1
Immunogen:	Fusion protein of human YY1AP1
Full name:	YY1 associated protein 1
Synonyms :	YAP; HCCA1; HCCA2; YY1AP
SwissProt:	Q9H869
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and Human brain

IHC Recommend dilution:

20-100

