

兔抗 ZBTB16 多克隆抗体

中文名称: 兔抗 ZBTB16 多克隆抗体

英文名称: Anti-ZBTB16 rabbit polyclonal antibody

别名: PLZF; ZNF145

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: ZBTB16

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene is a member of the Krueppel C2H2-type zinc-finger protein family and encodes a zinc finger transcription factor that contains nine Kruppel-type zinc finger domains at the carboxyl terminus. This protein is located in the nucleus, is involved in cell cycle progression, and interacts with a histone deacetylase. Specific instances of aberrant gene rearrangement at this locus have been associated with acute promyelocytic leukemia (APL). Alternate transcriptional splice variants have been characterized.
Applications:	ELISA, WB
Name of antibody:	ZBTB16

Immunogen:	Fusion protein of human ZBTB16
Full name:	zinc finger and BTB domain containing 16
Synonyms :	PLZF; ZNF145
SwissProt:	Q05516
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	74 kDa
WB Positive control:	Mouse fat and spleen tissue
WB Recommended dilution:	500-2000

