

兔抗 ZBTB18 多克隆抗体

中文名称: 兔抗 ZBTB18 多克隆抗体

英文名称: Anti-ZBTB18 rabbit polyclonal antibody

别名: RP58; MRD22; TAZ-1; ZNF238; C2H2-171

相关类别: 一抗

储存: 冷冻 (-20℃)

抗原: ZBTB18

宿主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a C2H2-type zinc finger protein which acts as a transcriptional repressor of genes involved in neuronal development. The encoded protein recognizes a specific sequence motif and recruits components of chromatin to target genes. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	ZBTB18
Immunogen:	Fusion protein of human ZBTB18
Full name:	zinc finger and BTB domain containing 18
Synonyms :	RP58; MRD22; TAZ-1; ZNF238; C2H2-171
SwissProt:	Q99592

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human lung cancer
IHC Recommend dilution:	20-100



