

兔抗 ZBTB40 多克隆抗体

中文名称：兔抗 ZBTB40 多克隆抗体

英文名称：Anti-ZBTB40 rabbit polyclonal antibody

别名：ZNF923

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZBTB40

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	WB, IHC
Name of antibody:	ZBTB40
Immunogen:	Synthesized peptide derived from C-terminal of human ZBTB40.
Full name:	zinc finger and BTB domain containing 40
Synonyms :	ZNF923
SwissProt:	Q9NUA8
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	138 kDa
WB Positive control:	293 cells lysate

WB Recommended dilution: 500-3000



