

## 兔抗 ZC3H4 多克隆抗体

中文名称：兔抗 ZC3H4 多克隆抗体

英文名称：Anti-ZC3H4 rabbit polyclonal antibody

级别：C19orf7

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZC3H4

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of CCCH (C-x8-C-x5-C-x3-H type) zinc finger domain-containing proteins. These zinc finger domains, which coordinate zinc finger binding and are characterized by three cysteine residues and one histidine residue, are nucleic acid-binding. Other family members are known to function in post-transcriptional regulation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZC3H4
<b>Immunogen:</b>	Synthetic peptide of human ZC3H4
<b>Full name:</b>	zinc finger CCCH-type containing 4
<b>Synonyms:</b>	C19orf7

<b>SwissProt:</b>	Q9UPT8
<b>ELISA Recommended dilution</b> :	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100

