

## 兔抗 ZDHHC5 多克隆抗体

中文名称: 兔抗 ZDHHC5 多克隆抗体

英文名称: Anti-ZDHHC5 rabbit polyclonal antibody

别名: DHHHC5; ZNF375

储存: 冷冻 (-20℃)

抗原: ZDHHC5

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat


标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5. Also palmitoylates FLOT2.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ZDHHC5
<b>Immunogen:</b>	Synthetic peptide of human ZDHHC5
<b>Full name:</b>	zinc finger, DHHC-type containing 5
<b>Synonyms:</b>	DHHHC5; ZNF375
<b>SwissProt:</b>	Q9C0B5
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	78 kDa

<b>WB Positive control:</b>	Mouse kidney tissue lysate
<b>WB Recommended dilution:</b>	200-1000



**mlbio** 梅联生物  
Good elisakit producers