

兔抗 ZFAND2B 多克隆抗体

中文名称: 兔抗 ZFAND2B 多克隆抗体

英文名称: Anti-ZFAND2B rabbit polyclonal antibody

别名: AIRAPL

抗原: ZFAND2B

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein containing AN1-type zinc-fingers and ubiquitin-interacting motifs. The encoded protein likely associates with the proteasome to stimulate the degradation of toxic or misfolded proteins. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, WB
Name of antibody:	ZFAND2B
Immunogen:	Synthetic peptide of human ZFAND2B
Full name:	zinc finger, AN1-type domain 2B
Synonyms :	AIRAPL
SwissProt:	Q8WV99

ELISA Recommended dilution:	1000-2000
WB Predicted band size:	28 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	200-1000



