

兔抗 ZFAND3 多克隆抗体

中文名称: 兔抗 ZFAND3 多克隆抗体

英文名称: Anti-ZFAND3 rabbit polyclonal antibody

别名: TEX27

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: ZFAND3

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	ZFAND3 contains 1 A20-type zinc finger and 1 AN1-type zinc finger. The exact function of ZF AND3 remains unknown.
Applications:	ELISA, IHC
Name of antibody:	ZFAND3
Immunogen:	Full length fusion protein
Full name:	zinc finger, AN1-type domain 3
Synonyms :	TEX27
SwissProt:	Q9H8U3
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

