

兔抗 ZFP2 多克隆抗体

中文名称: 兔抗 ZFP2 多克隆抗体

英文名称: Anti-ZFP2 rabbit polyclonal antibody

别名: ZNF751

抗原: ZFP2

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Probable transcription factor involved in neuronal differentiation and/or phenotypic maintenance.
Applications:	ELISA, IHC
Name of antibody:	ZFP2
Immunogen:	Fusion protein of human ZFP2
Full name:	ZFP2 zinc finger protein
Synonyms :	ZNF751
SwissProt:	Q6ZN57
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	30-150

