

## 兔抗 ZKSCAN3 多克隆抗体

中文名称：兔抗 ZKSCAN3 多克隆抗体

英文名称：Anti-ZKSCAN3 rabbit polyclonal antibody

别名：ZF47; Zfp47; ZFP306; ZNF306; ZNF309; zfp-47; ZSCAN13; ZSCAN35; dJ874C20.1; dJ874C20.1.

相关类别：一抗

抗原：ZKSCAN3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugat

### 技术规格

<b>Background:</b>	Transcriptional factor that binds to the consensus sequence 5'-[GT][AG][AGT]GGGG-3' and acts as a repressor of autophagy. Specifically represses expression of genes involved in autophagy and lysosome biogenesis/function such as MAP1LC3B, ULK1 or WIPI2. Associates with chromatin at the ITGB4 and VEGF promoters. Also acts as a transcription activator and promotes cancer cell progression and/or migration in various tumors and myelomas.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZKSCAN3

<b>Immunogen:</b>	Fusion protein of human ZKSCAN3
<b>Full name:</b>	zinc finger with KRAB and SCAN domains 3
<b>Synonyms :</b>	ZF47; Zfp47; ZFP306; ZNF306; ZNF309; zfp-47; ZSCAN13; ZSCAN35; dJ874C20.1; dJ874C20.1.
<b>SwissProt:</b>	Q9BRR0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	25-100

