

兔抗 ZNF8 多克隆抗体

中文名称：兔抗 ZNF8 多克隆抗体

英文名称：Anti-ZNF8 rabbit polyclonal antibody

别名：HF.18; Zfp128

抗原：ZNF8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	ZNF8 (zinc finger protein 8) is a 575 amino acid protein that belongs to the Krüppel C2H2-type zinc-finger protein family and contains seven C2H2-type zinc fingers and a KRAB domain. Localizing to the nucleus and ubiquitously expressed, ZNF8 is thought to play a role in transcriptional regulation and is encoded by a gene that maps to human chromosome 19q13.43.
Applications:	ELISA, IHC
Name of antibody:	ZNF8
Immunogen:	Fusion protein of human ZNF8
Full name:	zinc finger protein 8
Synonyms :	HF.18; Zfp128
SwissProt:	P17098

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human tonsil and human cervical cancer
IHC Recommend dilution:	30-150

