

兔抗 ZNF397 多克隆抗体

中文名称: 兔抗 ZNF397 多克隆抗体

英文名称: Anti-ZNF397 rabbit polyclonal antibody

别名: zinc finger protein 397; ZNF47; ZSCAN15

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: ZNF397

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein with a N-terminal SCAN domain, and the longer isoform contains nine C2H2-type zinc finger repeats in the C-terminal domain. The protein localizes to centromeres during interphase and early prophase, and different isoforms can repress or activate transcription in transfection studies. Multiple transcript variants encoding different isoforms have been found for this gene. Additional variants have been described, but their biological validity has not been determined.
Applications:	ELISA, IHC
Name of antibody:	ZNF397
Immunogen:	Fusion protein of human ZNF397

Full name:	zinc finger protein 397
Synonyms:	ZNF47; ZSCAN15
SwissProt:	Q8NF99
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human cervical cancer
IHC Recommend dilution:	40-200

