

兔抗 ZNF468 多克隆抗体

中文名称：兔抗 ZNF468 多克隆抗体

英文名称：Anti-ZNF468 rabbit polyclonal antibody

别名：zinc finger protein 468;

相关类别：一抗

抗原：ZNF468

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF468
Immunogen:	Full length fusion protein
Full name:	zinc finger protein 468
SwissProt:	Q5VIY5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	30-150

