

兔抗 ZNF483 多克隆抗体

中文名称: 兔抗 ZNF483 多克隆抗体

英文名称: Anti-ZNF483 rabbit polyclonal antibody

别名: ZSCAN48; ZKSCAN16

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: ZNF483

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF483
Immunogen:	Full length fusion protein
Full name:	zinc finger protein 483
Synonyms:	ZSCAN48; ZKSCAN16
SwissProt:	Q8TF39
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

