

兔抗 ZNF688 多克隆抗体

中文名称: 兔抗 ZNF688 多克隆抗体

英文名称: Anti-ZNF688 rabbit polyclonal antibody

别名: zinc finger protein 688

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: ZNF688

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, WB, IHC
Name of antibody:	ZNF688
Immunogen:	Fusion protein of human ZNF688
Full name:	zinc finger protein 688
SwissProt:	P0C7X2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and Human prostate cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	31 kDa
WB Positive control:	A172 and 231 cell lysates
WB Recommended dilution:	500-2000



