

小鼠抗 BRK1 单克隆抗体

中文名称: 小鼠抗 BRK1 单克隆抗体

英文名称: Anti-BRK1 mouse monoclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Mouse

抗 原: BRK1

反应种属: Human

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	Involved in regulation of actin and microtubule organization. Part of a WAVE complex that activates the Arp2/3 complex. As component of the WAVE1 complex, required for BDNF-NTRK2 endocytic trafficking and signaling from early endosomes.
Applications:	ELISA, WB
Name of antibody:	BRK1
Immunogen:	Fusion protein of human BRK1
Full name:	BRICK1, SCAR/WAVE actin nucleating complex subunit
Synonyms:	MDS027; hHBrk1; C3orf10; HSPC300
SwissProt:	Q8WUW1
WB Predicted band size:	9 kDa
WB Positive control:	H520 whole cell lysate
WB Recommended dilution:	200-1000

