

## 小鼠抗 CYCS 单克隆抗体

中文名称： 小鼠抗 CYCS 单克隆抗体

英文名称： Anti-CYCS mouse monoclonal antibody

别名： cytochrome c; somatic; CYC; HCS; THC4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： CYCS

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a small heme protein that functions as a central component of the electron transport chain in mitochondria. The encoded protein associates with the inner membrane of the mitochondrion where it accepts electrons from cytochrome b and transfers them to the cytochrome oxidase complex. This protein is also involved in initiation of apoptosis. Mutations in this gene are associated with autosomal dominant nonsyndromic thrombocytopenia. Numerous processed pseudogenes of this gene are found throughout the human genome.
<b>Applications:</b>	WB, IHC

<b>Name of antibody:</b>	CYCS
<b>Immunogen:</b>	Fusion protein of human CYCS
<b>Full name:</b>	cytochrome c, somatic
<b>Synonyms:</b>	CYC; HCS; THC4
<b>SwissProt:</b>	P99999
<b>IHC positive control:</b>	Human tonsil and Human thyroid cancer; Human ovarian cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	100-500
<b>WB Predicted band size:</b>	12 kDa
<b>WB Positive control:</b>	Hela, 231, HepG2 and Jurkat cell, Mouse kidney tissue lysates
<b>WB Recommended dilution:</b>	5000-20000



