

小鼠抗 ERP29 单克隆抗体

中文名称: 小鼠抗 ERP29 单克隆抗体

英文名称: Anti-ERP29 mouse monoclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Mouse

抗 原: ERP29

反应种属: Human

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a protein which localizes to the lumen of the endoplasmic reticulum (ER). It is a member of the protein disulfide isomerase (PDI) protein family but lacks an active thioredoxin motif, suggesting that this protein does not function as a disulfide isomerase. The canonical protein dimerizes and is thought to play a role in the processing of secretory proteins within the ER. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, WB
Name of antibody:	ERP29
Immunogen:	Fusion protein of human ERP29
Full name:	endoplasmic reticulum protein 29
Synonyms:	ERp28; ERp31; PDIA9; PDI-DB; C12orf8; HEL-S-107
SwissProt:	P30040

WB Predicted band size:	29 kDa
WB Positive control:	MCF-7 whole cell lysate
WB Recommended dilution:	1000-5000

