

## 兔抗 IL31RA 多克隆抗体

中文名称：兔抗 IL31RA 多克隆抗体

英文名称：Anti-IL31RA rabbit polyclonal antibody

别名：CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384

相关类别：一抗

储存：冷冻（-20℃）

抗原：IL31RA

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the type I cytokine receptor family. This receptor, with homology to gp130, is expressed on monocytes, and is involved in IL-31 signaling via activation of STAT-3 and STAT-5. It functions either as a monomer, or as part of a receptor complex with oncostatin M receptor (OSMR). Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IL31RA

<b>Immunogen:</b>	Fusion protein of human IL31RA
<b>Full name:</b>	interleukin 31 receptor A
<b>Synonyms :</b>	CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384
<b>SwissProt:</b>	Q8NI17
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human breast cancer and human prostate cancer
<b>IHC Recommend dilution:</b>	25-100

