

## 兔抗 KANSL1L 多克隆抗体

中文名称：兔抗 KANSL1L 多克隆抗体

英文名称：Anti-KANSL1L rabbit polyclonal antibody

别名：KAT8 regulatory NSL complex subunit 1 like; MSL1v2; C2orf67

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KANSL1L

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The function of this protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	KANSL1L
<b>Immunogen:</b>	Fusion protein of human KANSL1L
<b>Full name:</b>	KAT8 regulatory NSL complex subunit 1 like
<b>Synonyms:</b>	MSL1v2; C2orf67
<b>SwissProt:</b>	A0AUZ9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	100-300

