

## 兔抗 GTF3A 多克隆抗体

- 中文名称：兔抗 GTF3A 多克隆抗体
- 英文名称：Anti-GTF3A rabbit polyclonal antibody
- 别名：general transcription factor IIIA; AP2; TFIIIA
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：GTF3A
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

**Background:**

The product of this gene is a zinc finger protein with nine Cis[2]-His[2] zinc finger domains. It functions as an RNA polymerase III transcription factor to induce transcription of the 5S rRNA genes. The protein binds to a 50 bp internal promoter in the 5S genes called the internal control region (ICR), and nucleates formation of a stable preinitiation complex. This complex recruits the TFIIC and TFIIB transcription factors and RNA polymerase III to form the complete transcription complex. The protein is thought to be translated using a non-AUG translation initiation site in mammals based on sequence an

	alysis, protein homology, and the size of the purified protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GTF3A
<b>Immunogen:</b>	Synthetic peptide of human GTF3A
<b>Full name:</b>	general transcription factor IIIA
<b>Synonyms:</b>	AP2; TFIIIA
<b>SwissProt:</b>	Q92664
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	50-300

