

兔抗 GTF2H1 多克隆抗体

中文名称: 兔抗 GTF2H1 多克隆抗体

英文名称: Anti-GTF2H1 rabbit polyclonal antibody

级别: P62; BTF2; TFB1; TFIIH

储存: 冷冻 (-20℃)

抗原: GTF2H1

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | Component of the core-TFIIH basal transcription factor or involved in nucleotide excision repair (NER) of DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II. |
| Applications: | ELISA, IHC |
| Name of antibody: | GTF2H1 |
| Immunogen: | Synthetic peptide of human GTF2H1 |
| Full name: | general transcription factor IIH subunit 1 |
| Synonyms: | P62; BTF2; TFB1; TFIIH |
| SwissProt: | P32780 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human thyroid cancer |

IHC Recommend dilution:

20-100



