

兔抗 QSOX1 多克隆抗体

中文名称： 兔抗 QSOX1 多克隆抗体

英文名称： Anti-QSOX1 rabbit polyclonal antibody

别名： Q6; QSCN6

储存： 冷冻（-20℃）

抗原： QSOX1

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains domains of thioredoxin and ERV1, members of two long-standing gene families. The gene expression is induced as fibroblasts begin to exit the proliferative cycle and enter quiescence, suggesting that this gene plays an important role in growth regulation. Two transcript variants encoding two different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	QSOX1
Immunogen:	Fusion protein of human QSOX1

Full name:	quiescin sulphydryl oxidase 1
Synonyms:	Q6; QSCN6
SwissProt:	O00391
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	83 kDa
WB Positive control:	Human placenta tissue lysate
WB Recommended dilution:	200-1000



