

## 兔抗 RAB11A 多克隆抗体

中文名称：兔抗 RAB11A 多克隆抗体

英文名称：Anti-RAB11A rabbit polyclonal antibody

别名：YL8

相关类别：一抗

抗原：RAB11A

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RAB11A
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	RAB11A, member RAS oncogene family
<b>Synonyms :</b>	YL8

<b>SwissProt:</b>	P62491
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	Human fetal brain tissue , 231 cell
<b>WB Recommended dilution:</b>	500-2000



