

兔抗 RABEPK 多克隆抗体

中文名称: 兔抗 RABEPK 多克隆抗体

英文名称: Anti-RABEPK rabbit polyclonal antibody

别名: Rab9 effector protein with kelch motifs; p40; RAB9P40; bA65N13.1

抗原: RABEPK

储存: 冷冻 (-20°C)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Rab9 effector required for endosome to trans-Golgi network (TGN) transport.
Applications:	ELISA, IHC
Name of antibody:	RABEPK
Immunogen:	Fusion protein of human RABEPK
Full name:	Rab9 effector protein with kelch motifs
Synonyms:	p40; RAB9P40; bA65N13.1
SwissProt:	Q7Z6M1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human thyroid cancer
IHC Recommend dilution:	50-300

