

## 兔抗 KRTAP11-1 多克隆抗体

中文名称:	兔抗 KRTAP11-1 多克隆抗体
英文名称:	Anti-KRTAP11-1 rabbit polyclonal antibody
别名:	HACL1; HACL-1; KAP11.1
抗原:	KRTAP11-1
储存:	冷冻 (-20℃)
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human
标记物:	Unconjugate
克隆类型:	Unconjugate

### 技术规格

<b>Background:</b>	In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	KRTAP11-1
<b>Immunogen:</b>	Synthetic peptide of human KRTAP11-1
<b>Full name:</b>	keratin associated protein 11-1

<b>Synonyms :</b>	HACL1; HACL-1; KAP11.1
<b>SwissProt:</b>	Q8IUC1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	20-100



