

## 兔抗 RBM15 多克隆抗体

中文名称: 兔抗 RBM15 多克隆抗体

英文名称: Anti-RBM15 rabbit polyclonal antibody

别名: RNA binding motif protein 15; OTT; OTT1; SPEN

抗原: RBM15

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Members of the SPEN (Split-end) family of proteins, including RBM15, have repressor function in several signaling pathways and may bind to RNA through interaction with spliceosome components (Hiriart et al., 2005 [PubMed 16129689]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RBM15
<b>Immunogen:</b>	Fusion protein of human RBM15
<b>Full name:</b>	RNA binding motif protein 15
<b>Synonyms:</b>	OTT; OTT1; SPEN
<b>SwissProt:</b>	Q96T37
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human esophagus cancer and Human tonsil
<b>IHC Recommend dilution:</b>	30-150

