

## 兔抗 RBX1 多克隆抗体

中文名称：兔抗 RBX1 多克隆抗体

英文名称：Anti-RBX1 rabbit polyclonal antibody

别名：ROC1; RNF75; BA554C12.1

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：RBX1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This locus encodes a RING finger-like domain-containing protein. The encoded protein interacts with cullin proteins and likely plays a role in ubiquitination processes necessary for cell cycle progression. This protein may also affect protein turnover. Related pseudogenes exist on chromosomes 2 and 5.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	RBX1
<b>Immunogen:</b>	Fusion protein of human RBX1
<b>Full name:</b>	ring-box 1, E3 ubiquitin protein ligase
<b>Synonyms :</b>	ROC1; RNF75; BA554C12.1
<b>SwissProt:</b>	P62877

<b>IHC positive control:</b>	Human esophagus cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	12 kDa
<b>WB Positive control:</b>	Jurkat, MCF7, hela and hepg2 cell, mouse brain and human placenta tissue
<b>WB Recommended dilution:</b>	500-2000

