

## 兔抗 RCAN3 多克隆抗体

中文名称: 兔抗 RCAN3 多克隆抗体

英文名称: Anti-RCAN3 rabbit polyclonal antibody

别名: RCN3; MCIP3; hRCN3; DSCR1L2

抗原: RCAN3

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	RCAN3 inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. It could play a role during central nervous system development.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	RCAN3
<b>Immunogen:</b>	Synthetic peptide of human RCAN3
<b>Full name:</b>	RCAN family member 3
<b>Synonyms:</b>	RCN3; MCIP3; hRCN3; DSCR1L2
<b>SwissProt:</b>	Q9UKA8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	27 kDa

<b>WB Positive control:</b>	Mouse brain tissue lysate
<b>WB Recommended dilution:</b>	200-1000

