

兔抗 RCOR2 多克隆抗体

中文名称：兔抗 RCOR2 多克隆抗体

英文名称：Anti-RCOR2 rabbit polyclonal antibody

别名：REST corepressor 2;

相关类别：一抗

抗原：RCOR2

储存：冷冻（-20℃）

宿主：Rabbit



反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May act as a component of a corepressor complex that represses transcription.
Applications:	ELISA, IHC
Name of antibody:	RCOR2
Immunogen:	Synthetic peptide of human RCOR2
Full name:	REST corepressor 2
SwissProt:	Q8IZ40
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100



mlbio 梅联生物
Good elisakit producers